

# 2025 Quarterly Promotions



**Rodin RapidCeram & Glaze N2-Free**  
See Page 1



**Rodin Luting Cement**  
See Page 3



**SDF Pro**  
38% Silver Diamine Fluoride  
See Page 14



**ProMin F & ProMin**  
See Page 6



**AntiSplatr Flex Slim**  
See Page 7



**EzAim ComforTab**  
Adhesive Sensor Holder  
See Page 10





**RODIN**<sup>®</sup>  
3D Resin


A Pac-Dent Brand

# RAPIDCERAM<sup>™</sup>

Uncompromised Strength, Accelerated Speed



**NEW**  
— 2025 —



## RapidCeram

High-Speed Ceramic Nanohybrid Resin

Delivers up to 2X the speed of the Sculpture series

300 g/Bottle.....\$490.00

**Net Price: \$376.50**

OM3	<b>24057</b>
A1	<b>24058</b>
A2	<b>24059</b>
A3	<b>24060</b>



## Titan

Ceramic Nanohybrid Resin

Increased flexibility for long-lasting full arch restorations

1.2 kg/Bottle.....\$1,285.71

**Net Price: \$964.28**

B1	<b>23986</b>	A3	<b>23989</b>
A1	<b>23987</b>	OM1	<b>24003</b>
A2	<b>23988</b>	OM3	<b>23985</b>
C2	<b>23990</b>		



## Glaze N<sub>2</sub>-Free

Universal Light-Curing Sealant

Enables complete polymerization without the need of nitrogen curing

24031.....5 mL/Bottle

\$99.95 **Net Price: \$74.96**

24015.....15 mL/Bottle

\$249.99 **Net Price: \$187.49**



## Soft Reline

Silicone-Based Reline

Enables complete polymerization without the need of nitrogen curing



24015.....Soft Reline Kit

\$199.95 **Net Price: \$149.96**

1 x 10 ml VPS adhesive primer, 6 x wide mouth mixing tips, 1 x 50 ml VPS reline cartridge, 6 x applicator swabs, 1 x tapered cone finishing stone, 1 x brown polishing wheel

**Buy 1, Get 1 x 5 mL Bottle of Glaze N2-Free FREE!\***

**Buy 3, Get 1 FREE!\***  
(Mix & Match Different Shades)

**Buy 3, Get 1 FREE!\***

**Buy 3, Get 1 FREE!\***

## EnVision™ VENEER SYSTEM

**NEW**

- ✓ Streamlined digital workflow
- ✓ Best-in-class aesthetics
- ✓ Advanced ceramic-filled resin

### Two-Way Bond Strength

The Veneer Cement ensures a strong bilateral bond to the existing tooth structure as well as to the veneer.



- |              |                                     |                  |
|--------------|-------------------------------------|------------------|
| <b>24017</b> | Rodin EnVision Intro Kit- Shade 0M1 |                  |
| <b>24018</b> | Rodin EnVision Intro Kit- Shade 0M3 |                  |
| <b>24019</b> | Rodin EnVision Intro Kit- Shade A1  | (Reg. \$1469.00) |
| <b>24020</b> | Rodin EnVision Intro Kit- Shade B1  |                  |
| <b>24021</b> | Rodin EnVision Veneer Cement Kit    | (Reg. \$845.00)  |



**Buy 1, Get 1 x 5 mL bottle of Glaze N2-FREE!\***

**SAVE \$99.95**

## SCULPT SUCCESS WITH RODIN SCULPTURE 2.0

The **FIRST** and **ONLY** Crown and Bridge 3D Printed Resin with Over 60% Ceramic Filler Infused with Zirconia



**1.2 kg/Bottle** ..... \$1285.71      **Net Price: \$964.28**

- |              |     |              |    |              |    |
|--------------|-----|--------------|----|--------------|----|
| <b>23869</b> | 0M3 | <b>23875</b> | A1 | <b>23877</b> | A3 |
| <b>23871</b> | B1  | <b>23876</b> | A2 | <b>23881</b> | C2 |

**Buy 3, Get 1 FREE (Mix & Match)!\***



# LUTING CEMENT


Resin Modified Glass Ionomer Cement

## Proven Strength, Effortless Versatility.

**NEW!**



### Comparison of Bond Strength for Different Materials and Cements



\*\*Data provided by Dr. Russell Giordano of Boston University

### Best in Class Bond Strength on:

- ✓ 3D Resin based
- ✓ Glass-ceramic
- ✓ Zirconia Restorations

#### 451101 Rodin Luting Cement Starter Kit (Reg. \$599.99)

**Net Price: \$399.99**

- 2x 5g Clear Syringe
- 1x 5g White Syringe
- 1x 5g White Opaque Syringe
- 1x 5g Bleach Syringe
- 1x 3g Try-in Gel Clear Syringe
- 1x 3g Try-in Gel White Syringe
- 1x 3g Try-in Gel White Opaque Syringe
- 1x 3g Try-in Gel Bleach Syringe
- 1x 5 mL Rodin Bond Bottle (5ml)
- 2x PacEtch 1.5g Syringe
- 1x 5ml Silane Primer Bottle



#### Rodin Luting Cement Refill, 2 x 5g Syringes

- 451102 Clear (Reg. \$269.99)
  - 451103 White Syringes (Reg. \$269.99)
  - 451104 White Opaque (Reg. \$269.99)
- Net Price: \$179.99**

#### Rodin Try-In Gel Refill, 1 x 3g Syringes

- 451107 Clear (Reg. \$49.00)
- Net Price: \$32.67**

**2+1**

Buy 2, Get 1 of the Same FREE!\*


## RODIN BOND

Dental Adhesive System



### ALSO DELIVERS

- High Shear Bond Strength for conventional zirconia and composite restorations
- Maximum adhesive strength when coupled with SensiGuard desensitizer, ensuring patient comfort



Data from an Independent Research Institute\*\*

#### 471101 Rodin Bond Adhesive System Kit

- 1x Rodin Bond Adhesive 5mL
- 1x SensiGuard Desensitizer 5mL
- 2x Rodin ETCHPRO 1.5 gm Syringes
- 20x Dispensing tips for Rodin ETCHPRO
- 100x MicroPro Brushes
- 20x Mixing Wells
- 1x Technique Guide
- 1x Instruction for Use

(Reg. \$298.69) **NET PRICE: \$199.13**

#### 471103 Rodin Bond Unit Dose Kit

- 100x 0.1 mL Unit Doses

(Reg. \$360.49) **NET PRICE: \$240.33**

#### 481102 Rodin Bond System Unit Dose Kit

- 50x Rodin 0.1 mL Unit Doses
- 50x SensiGuard 0.1 mL Unit Doses

(Reg. \$315.17) **NET PRICE: \$210.12**

#### 471102 Rodin Bond Refill

- 5 mL bottle

(Reg. \$173.03) **NET PRICE: \$115.35**

#### 481103 SensiGuard Desensitizer

- 100x 0.1 mL Unit Doses

(Reg. \$279.12) **NET PRICE: \$186.08**

**2+1**

Buy 2, Get 1 of the Same FREE!\*

## IMATRIX™ CLEAR **NEW!**

### Sectional Matrix System

#### 411201 Starter Kit

- 1x Clear Green Standard Ring
- 1x Clear Yellow Narrow Ring
- 1x Polishing Kit
- 1x Forceps
- 1x Tweezers
- 35x Clear Yellow Small Wedges
- 35x Clear Blue Medium Wedges
- 35x Clear Orange Large Wedges
- 35x Clear Red 4.5 mm Matrices
- 35x Clear Blue 5.5 mm Matrices
- 35x Clear Green 6.5 mm Matrices

(Reg. \$496.53) **Net Price: \$372.40**

#### 411202 Intro Kit

- 1x Clear Green Standard Ring
- 1x Clear Yellow Narrow Ring
- 35x Clear Yellow Small Wedges
- 35x Clear Blue Medium Wedges
- 35x Clear Orange Large Wedges
- 35x Clear Red 4.5 mm Matrices
- 35x Clear Blue 5.5 mm Matrices
- 35x Clear Green 6.5 mm Matrices

(Reg. \$362.69) **Net Price: \$272.02**

#### 411204 2x Clear Green Standard Rings

#### 411205 2x Clear Yellow Narrow Rings

(Reg. \$204.29) **Net Price \$153.20**

#### 411206 Clear Red 4.5 mm Matrices, 100/pk

#### 411207 Clear Blue 5.5 mm Matrices, 100/pk

#### 411208 Clear Green 6.5 mm Matrices, 100/pk

(Reg. \$131.00) **Net Price: \$98.25**

#### 411209 Clear Yellow Small Wedge Refill, 100/pk

#### 411210 Clear Blue Medium Wedge Refill, 100/pk

(Reg. \$58.45) **Net Price: \$43.84**



#### 411212 1x Composite Polishing kit

(Reg. \$105.99) **Net Price: \$79.49**

#### PM-F1 1x iMatrix Forceps

#### PM-T1 1x iMatrix Pin Tweezers

(Reg. \$65.84) **Net Price: \$49.38**

### 3+1

Buy 3, Get 1 of the same FREE\*

## IMATRIX™

### Sectional Matrix System

#### PM-01 Starter Kit

- 1x Black Standard Ring
- 1x Gray Narrow Ring
- 1x Forceps
- 1x Tweezers
- 35x Yellow Small Wedges
- 35x Blue Medium Wedges
- 35x Orange Large Wedges
- 35x 3.5 mm Matrices
- 35x 4.5 mm Matrices
- 35x 5.5 mm Matrices
- 35x 6.5 mm Matrices

(Reg. \$355.16) **Net Price: \$266.37**

#### PM-02 Intro Kit

- 1x Black Standard Ring
- 1x Gray Narrow Ring
- 35x Yellow Small Wedges
- 35x Blue Medium Wedges
- 35x Orange Large Wedges
- 35x 3.5 mm Matrices
- 35x 4.5 mm Matrices
- 35x 5.5 mm Matrices
- 35x 6.5 mm Matrices

(Reg. \$257.45) **Net Price: \$193.09**

#### PM-R2 2x Black Standard Rings

#### PM-R3 2x Grey Narrow Rings

(Reg. \$142.90) **Net Price: \$107.18**

#### PM-M3 3.5 mm Matrices, 100/pk

#### PM-M4 4.5 mm Matrices, 100/pk

(Reg. \$54.20) **Net Price: \$40.65**

#### PM-WS Yellow Small Wedge Refill, 100/pk

#### PM-WM Blue Medium Wedge Refill, 100/pk

(Reg. \$29.07) **Net Price: \$21.80**



## IMATRIX™-DISPOSABLE **NEW!**

### Disposable Matrix Bands

Available in Mylar Band



Universal Shape, Stainless Steel, 50/pk

#### 411101 6 mm Wide, Blue

#### 411102 4.5 mm Narrow, Teal

Curved Shape, Stainless Steel, 50/pk

#### 411103 6 mm Wide, Orange

#### 411104 4.5 mm Narrow, Yellow

Universal Shape, Mylar, 50/pk

#### 411105 6 mm Wide, Red

(Reg. \$72.05) **Net Price: \$54.04**

### 3+1

Buy 3, Get 1 of the same FREE\*

### 3+1

Buy 3, Get 1 of the same FREE\*

## IMATRIX™ COATED **NEW!**

Coated Matrix Bands



- 411301** Dark Blue Matrices Refill, 3.5mm, 100/pk
- 411302** Red Matrices Refill, 4.5mm, 100/pk
- 411303** Blue Matrices Refill, 5.5mm, 100/pk
- 411304** Green Matrices Refill, 6.5mm, 100/pk
- 411305** Dark Gray Matrices Refill, 7.5mm, 100/pk

(Reg. \$131.00) **Net Price: \$98.25**

**3+1**

Buy 3, Get 1 of the same FREE!\*

## ETCHPRO™

38% Phosphoric Acid Etching Gel

One of the best products in 2024  
as rated by a non-profit  
research organization



- 135** 4x 1.2 mL Syringes  
(Reg. \$12.46) **Net Price : \$9.35**

**3+1**

Buy 3, Get 1 of the same FREE!\*

## OPTIFLOW™ II

Nano Flowable Composite



4x 1.5 gm Syringes with 20x Pre-bent needle Tips

- |                        |                      |
|------------------------|----------------------|
| <b>8201</b> Shade A1   | <b>8206</b> Shade B1 |
| <b>8202</b> Shade A2   | <b>8209</b> Shade B3 |
| <b>8203</b> Shade A3   | <b>8212</b> Shade C2 |
| <b>8204</b> Shade A3.5 | <b>8215</b> Shade D2 |

(Reg. \$38.06) **Net Price: \$30.45**

**4+1**

Buy 4, Get 1 of the same FREE!\*

## C&B PROMPT™ 10:1

Temporary Crown & Bridge Material



1x 50 mL Cartridge, 10x Tips

- PD123** Shade A1
- PD124** Shade A2
- PD125** Shade A3

(Reg. \$114.84) **Net Price: \$86.13**

**3+1**

Buy 3, Get 1 of the same FREE!\*

## PACTEMP™

Automix NE with Releasing Agent



- PD130**  
1x 6g Automix Syringe  
15x Brown Mixing Tips

(Reg. \$29.95) **Net Price: \$23.96**

**4+1**

Buy 4, Get 1 of the same FREE!\*

## PACSEAL™

Pit & Fissure Sealant LC



- PD-138**  
4x 1.2cc Syringes,  
Natural Shade 20x Brush Tips

(Reg. \$46.13) **Net Price: \$34.60**

**3+1**

Buy 3, Get 1 FREE!\*

# Experience Tomorrow's Dental Solution, Today.



## ProMin™ F & ProMin™

Tooth Desensitizing Paste with 0.20% Sodium Fluoride

Recommend a New Solution to Your Patients that:


- ✓ Can be applied with a toothbrush or custom tray
- ✓ Is non-invasive and easy for all ages

- ✓ Actively rebuilds enamel to relieve sensitivity with CPIC
- ✓ Prevents caries and white spot lesions
- ✓ Increase resistance to acidic challenges

## What is the CPIC technology?

**CPICs (calcium phosphate ion clusters)** are nanoscale clusters of calcium and phosphate ions that are specifically designed to penetrate the microscopic pores and defects in demineralized enamel, making them highly effective in repairing and strengthening teeth.

CPICs are engineered to mimic the body's natural mineralization processes, providing a more bioavailable and efficient source of minerals compared to traditional technologies.



No Fluoride

**341101**

ProMin Tooth Desensitizing Paste, 2 x 25 g tubes, 10 boxes/case  
(Reg. \$189.99) **Net Price: \$142.49**



with Fluoride

**341102**

ProMin F Tooth Desensitizing Paste with Sodium Fluoride, 2 x 25 g tubes, 10 boxes/case  
(Reg. \$189.99) **Net Price: \$142.49**



with Fluoride

**341103**

ProMin F Tooth Desensitizing Paste with Sodium Fluoride, 4.5 mL pen  
(Reg. \$19.99) **Net Price: \$14.99**

**3+1**

Buy 3, Get 1 of the Same FREE!\*

Patent Pending



# NEW! AntiSplatr® FLEX Slim

Disposable Prophy Angle

## Smart Flex Technology Adjusts to Every Tooth Shape



### Flex Cup

Easily adapts to the contours of any tooth



### Splatter Fender

Removes accumulated paste and bodily fluids



### Optimal Cup Lip

Greater flare and subgingival access



### 42% Slimmer Head

Better visibility and access

## AntiSplatr® Workflow

### AntiSplatr® Flex Slim

Disposable Prophy Angles

200/pk, \$169.95 **Net: \$127.46**

**212101** AntiSplatr Flex Slim

**212103** AntiSplatr Flex Ergo Slim

1200/pk, \$999.95 **Net: \$749.96**

**212102** AntiSplatr Flex Slim

NET PRICE  
as low as  
**62¢**  
per angle



Buy 3, Get 1 FREE! Or Buy 2 cases of 1200/case, Get 1 ProMate Contra Handpiece and 1 box of ProPaste One (Assorted, Coarse) FREE!\*

### ProMate™ Contra

Hygiene Contra Handpiece

**2** YEARS WARRANTY



Ergonomic Design

**251201**

ProMate Contra Hygiene Handpiece

\$799.00

**Net: \$599.25**



Buy 3, Get 1 FREE!\*

### ProPaste™ One

Disposable Prophy Ring & AntiSplatr Paste Combo

144/pk, \$65.91

**23101** Assorted, Coarse

**23102** Assorted, Medium

**23103** Mint, Coarse



**Net Price: \$43.94**

Buy 2, Get 1 of the Same FREE!\*



Scan to Request Your  
**FREE Samples Today**  
or visit [bit.ly/3z9R7i5](http://bit.ly/3z9R7i5)

\*Case Sensitive

## ANTISPLATR® Disposable Prophy Angle

- ✓ Compact Splatter Fenders effectively removes accumulated paste and bodily fluids
- ✓ Optimal Cup Lip Allows greater flare & subgingival access



100/pk, \$78.66

**LSFAS-100** Pink Soft  
**LSFAF-100** Blue Firm **Net Price: \$59.00**

500/pk, \$392.35

**LSFAS-500** Pink Soft  
**LSFAF-500** Blue Firm **Net Price: \$294.26**

**3+1**

Buy 3, Get 1 of the Same FREE!\*

## OSCILLATING ANTISPLATR® MINI Mini Disposable Prophy Angle



Oscillating motion effectively reduces splatter 20% slimmer body

NET PRICE  
as low as  
**70¢**  
per angle

**100/pk | ASMS-100**  
Soft Cup, Pink

(Reg. \$105.04) **Net Price: \$70.03**

**100/pk | ASMF-100**  
Firm Cup, White

(Reg. \$105.04) **Net Price: \$70.03**

**ASMB-100** Tapered Brush

(Reg. \$115.55) **Net Price: \$77.03**

**2+1**

Buy 2, Get 1 of the same FREE!\*

## PROMATE™ CL Cordless Hygiene Handpiece

ONE OF THE BEST CORDLESS  
PROPHY HANDPIECES OF 2023  
AS RATED BY A NON-PROFIT DENTAL  
RESEARCH ORGANIZATION



**PMCL-100**.....1 Handpiece Kit  
(Reg. \$1033.57) **Net Price: \$775.18**

## PROMATE™ EZ-Q Quick Disconnect Hygiene Handpiece



**PA-EZ-CB**.....Handpiece & Coupler  
(Reg. \$821.94) **Net Price: \$616.46**

## PROMATE™ LF Lube-Free Hygiene Handpiece

**PA-DM5H**.....1 Unit  
(Reg. \$607.06) **Net Price: \$455.30**

## PROMATE™ Hygiene Handpiece

**HH-100**.....1 Unit  
(Reg. \$567.52) **Net Price: \$425.64**

**3+1**

Buy 3 handpieces, Get 1 of the same FREE!\*

## FUNIMAL™ PEDO Pediatric Mini Disposable Prophy Angle

**100/pk**  
**167-100P** Assorted, Pedo Cup  
**167-100R** Assorted, Regular Cup  
(Reg. \$104.78) **Net Price: \$69.85**

NET PRICE  
as low as  
**70¢**  
per angle



**2+1**

Buy 2, Get 1 of the same FREE!\*

## MINI 90™/MINI ERGO™ Mini Disposable Prophy Angle

ONE OF THE BEST PRODUCTS OF 2018  
as rated by a non-profit dental research organization

**500/BX**  
**S163-90-500** 90° Soft Cup  
**S163-500** Ergo Soft Cup  
**S164-90-500** 90° Firm Cup  
**S164-500** Ergo Firm Cup  
(Reg. \$443.43) **Net Price: \$295.62**

NET PRICE  
as low as  
**59¢**  
per angle

 4.5 out of 5



**2+1**

Buy 2, Get 1 of the same FREE!\*

## PROANGLE® PLUS Disposable Prophy Angle

**144/BX**  
**249-144** Soft Teal Cup  
**249S-144** Torque Lavender Cup (Reg. \$71.34)  
**250-144** Firm Lavender Cup **Net Price: \$53.51**  
**251T-144** Super Soft Cup  
**100/BX**  
**238-100** Tapered Brush (Reg. \$84.00) **Net Price: \$63.00**



**3+1**

Buy 3, Get 1 of the same FREE!\*

# AUTOSCALER®

 **South East Instruments**  
A Pac-Dent Company



*Best-in-class warranty with 6 years for circuit boards and 2 years for all other parts!*

## AUTOSCALER® Ultrasonic Scaler

- ✓ 25K and 30K models available in tabletop or built-in configurations
- ✓ Made for special periodontal and any supra- subgingival or routine scaling
- ✓ Compatible with most industry-leading magnetostrictive insert tips
- ✓ Best-in-class warranty with 6 years for circuit boards and 2 years for all other parts
- ✓ Detachable handpiece assembly for maximum convenience
- ✓ Automatic tuning circuit provides even power under all load conditions
- ✓ Linear power control; allows the operator to choose the exact power level needed for all procedures




**AW25C** Tabletop Model, 25K.....\$1029.09  
**AW30C** Tabletop Model, 30K.....\$1029.09




**AW25BI** Build-In Model, 25K.....\$1655.13  
**AW30BI** Build-In Model, 30K.....\$1655.13

## AUTOSCALER® INSERT TIPS Ultrasonic Scaler Insert Tips

- ✓ Non-slip grip for greater tactile sensitivity
- ✓ The inserts use a base flow water delivery
- ✓ Compatible with most magnetostrictive handpieces
- ✓ 90 Day Warranty against defects arising from faulty material or workmanship

	<b>Universal</b>	Straight shank with medium length and tapered tip.
	<b>Deposit Types</b>	<ul style="list-style-type: none"> <li>• Light to heavy deposits</li> <li>• Tenacious calculus</li> </ul>
	<b>Applicable Areas</b>	<ul style="list-style-type: none"> <li>• Supragingival</li> <li>• Subgingival if tissue allows</li> </ul>
	<b>Power Setting</b>	Low to High

**INS 508601** 25k Ultrasonic Universal Insert Tip (TFI-10).....\$131.27  
**INS 508701** 30k Ultrasonic Universal Insert Tip (TFI-10).....\$131.27

	<b>Thin</b>	Straight shank with long length and a thin tapered tip.
	<b>Deposit Types</b>	Light to moderate deposits
	<b>Applicable Areas</b>	Subgingivally for periodontal pockets, but also can be used for supragingival applications
	<b>Power Setting</b>	Low to Medium

**INS 508602** 25k Ultrasonic Thin Insert Tip.....\$183.78  
**INS 508702** 30k Ultrasonic Thin Insert Tip.....\$183.78

## EZAIM® COMFORTAB™ **NEW!** Adhesive Cushioned Sensor Tab Holders



Compared to Dentsply Xios Plus and Schick AimRight™

- ✓ Cushioned tab for extra patient comfort
- ✓ Individually cut and packaged
- ✓ Strong adhesive allows for reliable positioning

**131121** Anterior Blue, 50/pk  
**131122** Anterior White, 50/pk  
**131123** Bitewing Red, 50/pk  
**131124** Bitewing Blue, 50/pk  
**131125** Yellow Posterior, 50/pk  
 (Reg. \$39.99) **Net Price: \$29.99**

**131115** Red SoftBite Tab, 100/pk  
 (Reg. \$50.99) **Net Price: \$38.24**

## Universal Adhesive Fit for Every Sensor



① Anterior



② Bitewing



③ Posterior



④ Cushioned Tab

**3+1**

Buy 3, Get 1 of the Same FREE!\*

## EZAIM® DEFEND 3-IN-1 **NEW!** X-Ray Collimator Positioning System



- ✓ Built-in collimator reduces radiation and improves imaging
- ✓ 3-in-1 color ring for bitewing, anterior, and posterior
- ✓ Compact design removes need for extra parts
- ✓ Compatible with all EzAim bite blocks for digital sensors, phosphor plates, and film
- ✓ Autoclavable

### EzAim Defend 3-in-1 Positioning System

**131127** Positioning System, 1/set  
 (Reg. \$149.50) **Net Price: \$112.13**

**131128** Positioning System, 2/set  
 (Reg. \$289.99) **Net Price: \$217.49**

**131129** Positioning Ring, 1/pk  
 (Reg. \$59.99) **Net Price: \$44.99**

**131102** Red Bitewing Positioning Arm, 1/pk  
 (Reg. \$44.91) **Net Price: \$33.68**

**131126** Blue/Yellow Positioning Arm, 1/pk  
 (Reg. \$44.91) **Net Price: \$33.68**

**131130** EzAim ComforTab Complete Intro Kit  
 1x Defend 3-in-1 Positioning Ring and 2x Positioning Arms,  
 60 ComforTab holders (10 of each: anterior, posterior,  
 bitewing, and bitetab)  
 (Reg. \$179.95) **Net Price: \$134.96**



①



②



③

- ① Anterior
- ② Bitewing
- ③ Posterior

**3+1**

Buy 3, Get 1 of the Same FREE!\*




## EZAIM® X-RAY POSITIONING SYSTEM

Compared to Dentsply Rinn XCP and Xios Plus

### Aiming Rings

Autoclavable



Blue Anterior

Red Bitewing



Yellow Posterior

Green Endo

- 131101** Blue Anterior Positioning Arm
- 131102** Red Bitewing Positioning Arm
- 131103** Yellow Posterior Positioning Arm
- 131104** Green Endo Positioning Arm  
(Reg. \$44.91)

- 131105** Blue Anterior Positioning Ring
- 131106** Red Bitewing Positioning Ring
- 131107** Yellow Posterior Positioning Ring
- 131108** Green Endo Positioning Ring  
(Reg. \$12.69)

Adhesive Sensor Holder

- 131109** Blue Anterior, 100/pk
- 131110** Yellow Posterior, 100/pk
- 131111** Red Bitewing, 100/pk
- 131112** GreenEndo Tab, 100/pk  
(Reg. \$70.02)

- 131114** EzAim X-Ray Positioning System Intro Kit  
(Reg. \$297.65)


Individually Cut Adhesive Sensor Holder **NEW!**

- 131116** Anterior Blue, 50/pk (Reg. 25.00)
- 131117** Anterior White, 50/pk
- 131118** Bitewing Red, 50/pk **Net Price: \$18.75**
- 131119** Bitewing Blue, 50/pk
- 131120** Yellow Posterior, 50/pk

### Positioning Arms

Autoclavable

- Compatible with existing aiming rings
- Autoclavable steel components



1	Yellow Posterior
2	Blue Anterior
3	Red Bitewing
4	Green Endo

### Sensor Holder Tabs

Disposable

- Peel-and-stick adhesive
- Universal Digital sensor compatibility

1	Red/Blue Bitewing
2	Yellow Posterior
3	Blue/White Anterior
4	Green Endo



**3+1**


Buy 3 Individually Cut Adhesive Sensor Tabs, Get 1 of the same FREE!\*

## SOFTEDGE™ CUSHIONS **NEW!**

X-Ray Film and Sensor Cushions

Compare to Crosstex® Edge-Ease® Cushions

- 131301** Large (Blue) Fits Long Edge of Size 2 Film and Sensors, 300/bx
- 131302** Small (White) Fits Short Edge of Size 1-2 Film Size and All 0-1 Sensors, 300/bx  
(Reg. \$51.49) **Net Price: \$38.62**



Easy Peel



**3+1**

Buy 3 of the same, Get 1 FREE!\*

## EASYPROP™ **NEW!**

Disposable Mouth Props

Compared to EZ-Prop®, Stabe® & Supa Film Holders, Dentsply Rinn XCP® Bite Block

- 131201** EasyProp™ Disposable Mouth Props, 100 pcs/box
- 131202** EasyProp™ Disposable X-Ray Film Holder, 100 pcs/box
- 131203** EasyProp™ Disposable Bite Blocks, 100 pcs/box  
(Reg. \$59.33) **Net Price: \$44.50**



**3+1**

Buy 3 of the same, Get 1 FREE!\*

## SURETIGHT™

Air/Water Syringe Tips  
Comparable to Kerr Seal-Tight®

**109-200W** SureTight White/w Yellow Core, 200/pk  
**109-200AST** SureTight Colors  
**541101** SureTight Clear  
 (Reg. \$38.57) **Net Price: \$28.93**

**109-1500W** SureTight White/w Yellow Core, 1500/pk  
**109-1500AST** SureTight Colors  
**541102** SureTight Clear  
 (Reg. \$267.86) **Net Price: \$200.90**



**3+1**

Buy 3, Get 1 of the Same FREE!\*

## NEOTIP™

Air/Water Syringe Tips  
Comparable to Flash Tips® & Crystal Tips®

**105-250AST** Assorted Colors, 250/pk  
 (Reg. \$52.40) **Net Price: \$39.30**

**105-1500AST** Assorted Colors, 1500/pk  
 (Reg. \$285.65) **Net Price: \$214.24**



**3+1**

Buy 3, Get 1 of the Same FREE!\*

## TRUTIP™ PLUS COLORS

Air/Water Syringe Tips

**108-150C** Clear, 150/pk  
 (Reg. \$38.96) **Net Price: \$29.22**

**108-150AST** Assorted Colors, 150/pk  
 (Reg. \$43.02) **Net Price: \$32.27**

**108-1500C** Clear, 1500/pk  
 (Reg. \$357.18) **Net Price: \$267.89**

**108-1500AST** Assorted Colors, 1500/pk  
 (Reg. \$384.15) **Net Price: \$288.11**



**3+1**

Buy 3, Get 1 of the Same FREE!\*

## Armored against Cross-Contamination with Time- and Space-Saving Solutions






## ARMOR™ COMBO

Air / Water Syringe Tip & Sleeve Combo

One of the best products in 2019 as rated by non-profit research organization



-  Saves time and space as a single item with single step
-  Better prevents cross-contamination with a tight fit and tactile-feel sleeve
-  With one-box packaging, reduces waste and saves valuable operating space



- 101-144** A/W Syringe Tips & Sleeve Combo, 144/pk  
**NEW!**  
**541201** A/W SureTight Syringe Tips & Sleeve Combo, 144/pk  
 (Reg. Price: \$56.77)      **Net Price: \$37.85**

**2+1**

Buy 2, Get 1 of the Same FREE!\*

# ENDODONTIC

## PACENDO™ Irrigation Solution

Each bottle comes with 1x luer-lock dispensing cap



Luer-Lock Bottle Cap

- |                  |                                    |                |                           |
|------------------|------------------------------------|----------------|---------------------------|
| <b>PE-ENC-16</b> | 2-in-1 EDTA & Chlorhexidine, 16 oz | (Reg. \$79.04) | <b>Net Price: \$59.28</b> |
| <b>PE-ENC-4</b>  | 2-in-1 EDTA & Chlorhexidine, 4 oz  | (Reg. \$28.66) | <b>Net Price: \$21.50</b> |
| <b>PE-SHB</b>    | 6% Sodium Hypochlorite, 16 oz      | (Reg. \$41.19) | <b>Net Price: \$30.89</b> |
| <b>PE-EDB</b>    | 17% EDTA, 16 oz                    | (Reg. \$42.85) | <b>Net Price: \$32.14</b> |
| <b>PE-ENC-16</b> | 95% Ethyl Alcohol, 4 oz            | (Reg. \$8.45)  | <b>Net Price: \$6.34</b>  |
| <b>PE-CHXB</b>   | 2% Chlorhexidine, 16 oz            | (Reg. \$44.60) | <b>Net Price: \$33.45</b> |



**3+1**

Buy 3, Get 1 of same item Free!\*

## PACENDO™ VITAL-ICE™ Pulp Vitality Spray

- PE-VI-01** Pulp Vitality Spray, 6 oz  
 (Reg. \$26.88)      **Net Price: \$20.16**



**3+1**

Buy 3, Get 1 FREE!\*

## SDF PRO™ **NEW!** 38% Silver Diamine Fluoride Solution

Compared to Advantage Arrest® Silver Diamine Fluoride


- ✓ Reduces dentinal hypersensitivity and prevents decay on high-risk sites
- ✓ Penetrates healthy dentin and enamel without discoloration
- ✓ Non-invasive—no drills, needles, or anesthetic needed
- ✓ Ideal when traditional treatments like fillings or sealants aren't an option



### How to Apply SDF PRO™




Remove excess saliva and isolate the area with cotton or gauze. Apply petroleum jelly if near the gingiva.



Thoroughly dry the affected tooth surface(s).



Dispense 1 drop of SDF Pro into a dappen dish and apply with an applicator.



Allow SDF Pro to absorb for up to one minute. Air dry to remove any excess saliva and precipitation.

**73014** 5 mL Bottle Kit: (Reg. \$119.95)  
1x 5 mL Bottle  
50 x Micro Applicator Tips **Net Price: \$89.99**

**73013** Unit Dose Kit: (Reg. \$139.99)  
30 x 0.1 mL Unit Doses  
50 x Micro Applicator Tips **Net Price: \$104.99**

**3+1**

Buy 3, Get 1 of the Same FREE!\*

## PERIO ARMOR® **NEW!** 1.7% HP Oral Cleansing Gel

ONE OF THE BEST PRODUCTS OF 2023  
AS RATED BY A NON-PROFIT DENTAL  
RESEARCH ORGANIZATION

- ✓ Features 1.7% hydrogen peroxide formula
- ✓ Non-invasive periodontal tray solution
- ✓ Helps protect against gum disease
- ✓ Helps eliminate harmful bacteria

**93101** Mint, 3.0 oz/tube, 6 tubes/pk  
**93102** Strawberry, 3.0 oz/tube, 6 tubes/pk  
**93103** Bubblegum, 3.0 oz/tube, 6 tubes/pk  
(Reg. \$93.68) **Net Price: \$70.26**



**3+1**

Buy 3, Get 1 of the Same FREE!\*

\*To qualify and receive your promotional goods, please submit your dealer invoice to [info@gingi-pak.com](mailto:info@gingi-pak.com). Cannot be combined with any other offer. No Cash Value. Equal or lesser value item(s) only. Promotions are subject to change or cancellation without notice. Please allow 8-10 weeks for delivery.

## GINGICAINE® IN SYRINGE

20% Benzocaine Topical Anesthetic

One of the best products in 2018  
as rated by a non-profit  
research organization



**20118S** Strawberry, 20x 1.2 mL Syringes  
(Reg. \$35.66) **Net Price: \$26.75**

**3+1**

Buy 3, Get 1 FREE!\*

## GINGICAINE® IN JAR

20% Benzocaine Topical Anesthetic

**20108** Bing Cherry, 1 oz jar  
**20114** Chocolate Mint, 1 oz jar  
(Reg. \$10.06) **Net Price: \$7.19**

**20109** Bing Cherry, 1 oz jar, 6/pk  
**20111** Assorted, 6/pk  
(Reg. \$58.09) **Net Price: \$41.49**



**5+2**

Buy 5, Get 2 of the Same FREE!\*

## Z TWIST® WEAVE

Braided Retraction Cords



**1117TM** Non-Impregnated Cord, #000 (Extremely Thin)  
**11170M** Non-Impregnated Cord, #00 (Very Thin)  
**10171M** Cord w/ Epinephrine, #1 (Thin)  
**10170M** Cord w/ Epinephrine, #00 (Very Thin)  
**12172M** Cord w/ Aluminum Sulfate, #2 (Medium)  
(Reg. Price: \$14.83 - 15.69) **Net Price: \$12.36 - \$13.08**

**5+1**

Buy 5, Get 1 of the Same FREE!\*

## GINGIPASTE™

15% Aluminum Chloride Retraction Paste



**12180** Intro Pack:  
24x Unit Doses  
1x Dispensing Gun (Reg. \$139.27) **Net Price: \$104.45**

**12181** Unit Dose Pack:  
24x Unit Doses (Reg. \$100.30) **Net Price: \$75.23**

**12182** 7x Syringe Pack:  
7x 0.7 gm Syringes  
15x Dispensing Tips (Reg. \$75.96) **Net Price: \$56.97**

**3+1**

Buy 3, Get 1 of the Same FREE!\*

## GINGI-AID® & STASIS®

Hemostatics



**12115** Gingi-Aid 25% Aluminum Chloride, 40 mL (Reg: \$46.22)  
**13115** Stasis 21% Ferric Sulfate, 40 mL **Net Price: \$34.67**  
**12116** Gingi-Aid Gel in Syringe, 1.2 mL, 4/pk (Reg: \$20.00)  
**13116** Stasis Gel in Syringe, 1.2 mL, 4/pk **Net Price: \$15.00**

**3+1**

Buy 3, Get 1 of the Same FREE!\*

## KROMATICA®

Color-Changing Alginate Impression Material

10 lbs/case  
**305-101** Fast Set, Cherry  
**305-102** Fast Set, Mint  
**305-103** Fast Set, Bubblegum  
**305-104** Fast Set, Unflavored  
**305-105** Normal Set, Cherry  
**305-106** Normal Set, Mint  
**305-107** Normal Set, Bubblegum



(Reg: \$143.41) **Net Price: \$119.51**



**5+1**

Buy 5, Get 1 of the Same FREE!\*

\*To qualify and receive your promotional goods, please submit your dealer invoice to [info@gingi-pak.com](mailto:info@gingi-pak.com). Cannot be combined with any other offer. No Cash Value. Equal or lesser value item(s) only. Promotions are subject to change or cancellation without notice. Please allow 8-10 weeks for delivery.